



MG-NA is a RNA/DNA extraction kit for isolation and purification of nucleic acids from tissues, whole blood, serum, plasma and nasopharyngeal swabs.

Kit Components	Volumes for 50 reaction	Storage conditions
Lysis Buffer (LB) (shake before use)	23 ml	Room Temperature
B-Mercaptoethanol (BME)	1 ml	Room Temperature
Binding Buffer (BB)	30 ml	Room Temperature
Magnetic Nano-beads (MNBs)	4 ml	Room Temperature
Washing Buffer I (WB1)	38 ml	Room Temperature
Washing Buffer II (WB2)	75 ml	Room Temperature
Elution Buffer	12 ml	Room Temperature

### Extraction Protocol

**(NB: For nucleic acid extraction from tissues, 200 mg tissue is homogenised in 1000  $\mu$ l lysis buffer, incubate at room temperature for 5 min, moving directly to step 3)**

1. Add 10  $\mu$ l BME to each 1 ml of lysis buffer (shake before use) directly before use.
2. Add 200  $\mu$ l of sample to 300  $\mu$ l lysis buffer in a 2 ml tube  $\rightarrow$  Tube 1. Mix by pipetting 3–5 times up and down. Incubate at room temperature for 5 minutes.
3. Shake and mix the magnetic beads well before use. Add 50  $\mu$ l of MNBs to Tube 1 the and mix well by pipetting up and down gently.
4. Add 400  $\mu$ l of binding buffer to Tube 1, then mix gently.
5. Incubate Tube 1 for 3 minutes at 65°C.
6. Insert Tube 1 into the magnetic rack and allow it to stand for 3 minutes at room temperature. Discard the clear solution, using a pipette.
7. Add 500  $\mu$ l washing buffer 1 to Tube 1 and repeat step 6.
8. Add 500  $\mu$ l washing buffer 2 to Tube 1 and repeat step 6. Repeat wash 2 step.
9. Remove Tube 1 from magnetic rack and resuspend MNBs in 50  $\mu$ l elution buffer (**Elute in 150  $\mu$ L in case of Genomic DNA tests**) and incubate it for 3 minutes at 65°C.
10. Apply magnetic separation and collect the RNA/DNA using a pipette while hot.

### Equipment needed

- 1- Magnetic Rack
- 2- Thermal Block

